


## Product Datasheet

### PTEN (Phospho-Ser370) antibody (orb106181)

<b>Description</b>	Rabbit polyclonal to PTEN (Phospho-Ser370). PTEN a phosphoinositide 3-phosphatase (PIP3) and a tumor...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB:1:500-1:2000, IHC-P:1:50-1:200
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HeLa whole cell...