


Product Datasheet

AFX (Phospho-Ser197) antibody (orb106196)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal to AFX (Phospho-Ser197). FOXO4 is a transcription factor AFX1 containing 1 fork-he... |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Dilution Range | WB:1:500-1:2000, IHC-P:1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of 293 whole cell ...