

Product Datasheet

G3BP-1 (Phospho-Ser232) antibody (orb106238)

Description Rabbit polyclonal to G3BP-1 (Phospho-Ser232). G3BP-1 an hnRNA-binding protein and endoribonuclease t...

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of 293 whole cell ...