



Product Datasheet

CD44 antibody (orb10628)


Description	Rabbit polyclonal antibody to CD44
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein of human CD44(NP_000601.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80KDa
Uniprot ID	P16070
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF/ICC: 1:50-1:100, FACS: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human tissue.