


Product Datasheet

SARS-CoV-2 N antibody (orb1063082)

Description	SARS-CoV-2 N antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Target	SARS-CoV-2 N
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	PODTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-SARS-CoV-2 N antibody