



Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized antibody (orb1063822)

| | |
|----------------------------|---|
| Description | SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized antibody |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant SARS-CoV-2 Nucleocapsid Protein (1-419/419 aa) |
| Target | SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized |
| Preservatives | 0.01M PBS(pH7.4) |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 46 kDa |
| Uniprot ID | P0DTC2 |
| Dilution Range | WB=1:500-2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis using Anti-SARS-CoV-2



Western blot analysis using Anti-SARS-CoV-2