

## Product Datasheet

# SARS-CoV-2 (2019-nCoV) Spike RBD (WT & Omicron) antibody (orb1063825)

<b>Description</b>	SARS-CoV-2 (2019-nCoV) Spike RBD (WT & Omicron) antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant SARS-CoV-2 Spike S1 Protein (14-685/1273 aa)
<b>Target</b>	SARS-CoV-2 (2019-nCoV) Spike RBD (WT & Omicron)
<b>Preservatives</b>	0.01M PBS(pH7.4)
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot Dependent
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>MW</b>	140 kDa
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Expiration Date</b>	12 months from date of receipt.