


Product Datasheet

Rat MCP-1 ELISA Kit (orb106599)

Description	Rat MCP-1 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q24-N148
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles. (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:5, concentration in serum is 300-900 pg/ml
Sample Types	cell culture supernatants and serum.
MW	16460 MW
Uniprot ID	P14844
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat MCP-1/CCL2 PicoKine ELISA Kit standa...