


Product Datasheet

Human Endothelin 1/EDN1 ELISA Kit (orb106603)

Description	Human Endothelin 1/EDN1 ELISA Kit
Reactivity	Human
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	A peptide reflecting the target protein in human. Sequence: C53-W73
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	24425 MW
Uniprot ID	P14138, P20800, P05305
Sensitivity	0.5 pg/ml
Expiration Date	Please enquire.



Human Endothelin PicoKine ELISA Kit stan...