


Product Datasheet

Mouse Angiopoietin-2 ELISA Kit (orb106607)

Description	Mouse Angiopoietin-2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: M1-F496
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 10-20 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	56576 MW
Uniprot ID	O35608
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse Angiopoietin-2 PicoKine ELISA Kit ...