


Product Datasheet

TNNI2 antibody (orb106870)

Description	Goat polyclonal antibody to TNNI2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TNNI2 (aa83-95)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 25kDa band observed in Human, Mouse and Rat Skeletal Muscle lysates (calculated MW of 21.3kDa according to NP_003273.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EVRVQKTSKELED
Clonality	Polyclonal
MW	21.3; 21.4
Entrez	7136
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human, Mouse an...