

Product Datasheet

Creatine Kinase BB antibody (orb106883)

Description	Mouse monoclonal to Creatine Kinase BB. Reversibly catalyzes the transfer of phosphate between ATP a...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	CKB
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	42644
Uniprot ID	P12277
NCBI	NP_001814.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

