

Product Datasheet

MCM9 antibody (orb106897)

Description	Rabbit polyclonal to MCM9. MCM9, a member of the MCM2-8 family, binds to chromatin and is required f...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MCM9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	127313
Uniprot ID	Q9NXL9
NCBI	NP_060166.2 , NP_694987.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

