

Product Datasheet

HHV14 UL38 antibody (orb106935)

Description	Rabbit polyclonal to HHV14 UL38. Structural component of the T=16 icosahedral capsid. The capsid is ...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TRX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50232
Uniprot ID	P17586
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

