

Product Datasheet

INSR antibody (orb106939)

Description	Rabbit polyclonal to INSR. Receptor tyrosine kinase which mediates the pleiotropic actions of insuli...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	INSR
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	156333
Uniprot ID	P06213
NCBI	NP_001073285.1 , NP_000199.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

