

Product Datasheet

FRA10AC1 antibody (orb106973)

Description	Rabbit polyclonal to FRA10AC1. The function of the C10orf4 protein is unknown. Expansion of a polymo...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FRA10AC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37548
Uniprot ID	Q70Z53
NCBI	NP_660289.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

