

Product Datasheet

L3MBTL3 antibody (orb106977)

Description	Rabbit polyclonal to L3MBTL3. Putative Polycomb group (PcG) protein. PcG proteins maintain the trans...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	L3MBTL3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88337
Uniprot ID	Q96JM7
NCBI	NP_115814.1 , NP_001007103.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

