

Product Datasheet

CUTA antibody (orb106980)

Description	Rabbit polyclonal to CUTA. May forms part of a complex of membrane proteins attached to acetylcholin...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CUTA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19116
Uniprot ID	O60888
NCBI	NP_057005.1 , NP_001014837.1 , NP_001014840.1 , NP_001014838.1 , NP_001014433.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

