

Product Datasheet

IFNA8 antibody (orb107004)

Description	Rabbit polyclonal to IFNA8. Produced by macrophages, IFN-alpha have antiviral activities. Interferon...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IFNA8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21989
Uniprot ID	P32881
NCBI	NP_002161.2
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.

