

Product Datasheet

BCL2 (Phospho-S70) antibody (orb107109)

Description	Rabbit polyclonal antibody to BCL2 (Phospho-S70)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	BCL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26266
Uniprot ID	P10415
NCBI	NP_000624.2
Dilution Range	DOT: 1:500, WB: 1:1000
Expiration Date	12 months from date of receipt.

