


Product Datasheet

NFκB-p105/p50(Phospho-Ser893) Antibody (orb107127)

Description	NFκB-p105/p50(Phospho-Ser893) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P19838
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...