


Product Datasheet

BCL-2(Phospho-Thr56) Antibody (orb107128)

Description	BCL-2(Phospho-Thr56) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P10415
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...