



Product Datasheet

c-Abl(Phospho-Tyr412) Antibody (orb107131)

Description	c-Abl(Phospho-Tyr412) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P00519
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Immunohistochemical analysis of formalin...