



## Product Datasheet

# AMPK $\alpha$ 1/AMPK $\alpha$ 2(Phospho-Thr183/Thr172) Antibody (orb107132)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | AMPK $\alpha$ 1/AMPK $\alpha$ 2(Phospho-Thr183/Thr172)<br>Antibody                        |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IHC, WB   |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Recombinant   |
| <b>Uniprot ID</b>          | <a href="#">Q13131</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Immunohistochemical analysis of formalin...



Western blot analysis of extracts from C...