



Product Datasheet

GRASP1 antibody (orb107143)

Description	GRASP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q4V328
Expiration Date	12 months from date of receipt.



Western blot analysis of A: HepG2 cells ...



Confocal immunofluorescence analysis of ...