


## Product Datasheet

### Mouse Endothelin 1/EDN1 ELISA Kit (orb107176)

<b>Description</b>	Mouse Endothelin 1/EDN1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	3.9 pg/ml - 250 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	A peptide reflecting the target protein in human. Sequence: C53-W73
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	22809 MW
<b>Uniprot ID</b>	<b>P22387</b>
<b>Sensitivity</b>	0.5 pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse Endothelin PicoKine ELISA Kit stan...