


Product Datasheet


CDKN1A Antibody (orb107328)

Description	CDKN1A Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P38936
Expiration Date	12 months from date of receipt.



Western blot analysis of Gel: 8%+10%SDS-

...



Immunohistochemical analysis of formalin...