



Product Datasheet

ANXA6 Antibody (orb107344)

Description	ANXA6 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08133
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of 10%SDS-PAGE lys...