



Product Datasheet

CBX8 Antibody (orb107362)

Description	CBX8 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9HC52
Expiration Date	12 months from date of receipt.



Western blot analysis of 231, HeLa, A549...



Immunohistochemical analysis of formalin...