



Product Datasheet

CLEC4D Antibody (orb107373)

| | |
|----------------------------|---|
| Description | CLEC4D Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | Q8WXI8 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Human fetal kid...



Immunohistochemical analysis of formalin...