



Product Datasheet

TEKT1 Antibody (orb107447)

Description	TEKT1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q969V4
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of NIH/3T3 cell ly...