



Product Datasheet

CUL4A Antibody (orb107489)

Description	CUL4A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13619
Expiration Date	12 months from date of receipt.



Western blot analysis of Human fetal mus...



Immunohistochemical analysis of formalin...