



Product Datasheet

GRO beta antibody (orb10749)


Description	GRO beta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MIP2 (51-107/107 aa)
Target	CXCL2/GRO Beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12 kDa
Uniprot ID	P19875
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=0.2ug/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Plasma Lysa...



Immunohistochemical staining of human re...



Immunohistochemical staining of mouse ut...