



Product Datasheet

F2R Antibody (orb107570)

Description	F2R Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P25116
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of K562 cel lysate...