



Product Datasheet

PDK4 Antibody (orb107572)

Description	PDK4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16654
Expiration Date	12 months from date of receipt.



Western blot analysis of Gel: 8%SDS-PAGE...



Immunohistochemical analysis of formalin...