
Product Datasheet

Cyclin-dependent kinase 7 Cdk7 Antibody (orb107628)

Description

Cyclin-dependent kinase 7 Cdk7 Antibody


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western blot analysis using Anti-Cdk7 an...

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of mouse Cdk7 (329-346aa RKRAEALEQGILPKKLIF), identical to the related rat sequence, and different from the related human sequence by two amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

38968 MW

Uniprot ID
[Q03147](https://www.uniprot.org/entry/Q03147)
Expiration Date

12 months from date of receipt.