

Product Datasheet

CD134/OX40/TNFRSF4 Antibody (orb107639)

Description

CD134/OX40/TNFRSF4 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of mouse CD134 (22-43aa ARRLNCVKHTYPSGHKCCRECQ).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for TNFRSF4 is approximately 0.5ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

30154 MW

Uniprot ID

P47741

Expiration Date

12 months from date of receipt.

1 2 3 — 250KD

— 100KD
— 70KD
— 55KD
— 35KD
— 25KD
— 15KD

Western blot analysis of
Lane 1:
Mouse B...