

Product Datasheet

CD134/OX40/TNFRSF4 Antibody (orb107640)

Description

CD134/OX40/TNFRSF4 Antibody

Species/Host

Rabbit

Reactivity

Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of rat CD134 (23-42aa KLNCVKDTYPSGHKCCRECQ).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for TNFRSF4 is approximately 5ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality


Polyclonal

MW

29895 MW

Uniprot ID
P15725
Expiration Date

12 months from date of receipt.



Western blot analysis of
Lane 1: Rat
Bra...