


Product Datasheet

Rat BD-3 protein (orb107689)

Description	Beta defensin-3, also known as BD-3 and DEFB-3, is a membrane active cationic peptide that functions...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	BD-3
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	KKVYNAVSCM TNGGICWLKC SGTFREIGSC GTRQLKCKK K
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	4.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat BD-3 protein