


Product Datasheet

Murine Betacellulin protein (orb107691)

Description	Recombinant Murine Betacellulin. Betacellulin (BTC) encoded by the BTC gene located on the chromosom...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Betacellulin
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	DGNTTRTPET NGS LCGAPGE NCTGTTPRQK VKTHFSRCPK QYKHYCIHGR CRFVDEQTPSCICEKGYFG ARCERVDLFY
Purity	> 96% by SDS-PAGE and HPLC analyses
MW	9.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine Betacellulin...