

Product Datasheet

Rat CINC-2 beta protein (orb107697)

| | |
|----------------------------|--|
| Description | The three GRO cDNAs encode 107 amino acid precursor proteins from which the N-terminal 34 amino acid... |
| Conjugation | Unconjugated |
| Tested Applications | HPLC, SDS-PAGE |
| Target | CINC-2 beta |
| Preservatives | Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, pH 7.4, 50mM NaCl |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only |
| Protein Sequence | RELRCQLKT LPRVDFENIQ SLTVTPPGPH CTQTEVIATL KDGQEVCLNP QAPRLQKIIQKLLKSDKSS |
| Purity | > 95% by SDS-PAGE and HPLC analyses |
| MW | 7.8 kDa |
| Expiration Date | 6 months from date of receipt. |