


Product Datasheet

Rat CNTF protein (orb107699)

Description	Recombinant Rat Ciliary Neurotrophic Factor. Ciliary neurotrophic factor (CNTF) is a polypeptide hor...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	CNTF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4, 2mM C-mercaptoethanol.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	AFAEQTPLTL HRRDLCSRSI WLARKIRSDL TALMESYVKH QGLNKNINLDSVDGVPVAST DRWSEMTEAE RLQENLQAYR TFQGMLTKLL EDQRVHFTPTGDFHQAIHT LMLQVSAFAY QLEELMVLLLE QKIPENEADG MPATVGDGGLFEKKLWGLKV LQELSQWTVR SIHDLRVISS HQMGISALES HYGAKDKQM
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	22.7 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat CNTF protein