


Product Datasheet

Cardiotrophin-1 protein (orb107701)

Description	Recombinant Rat CT-1 is a member of the IL-6 family of cytokines which also includes LIF, CNTF, OSM ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Cardiotrophin-1
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Cardiotrophin-
1 pro...