

Product Datasheet

Rat CXCL17 protein (orb107702)

Description	Recombinant Rat VEGF Co-regulated Chemokine 1/CXCL17. CXCL17 also named Dendritic Cell and Monocyte ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	CXCL17
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SPNQEVARHH GDQHQAPRRW LWEGGQECDC KDWSLRVSKR KTTAVLEPPRKQCPCDHVKG SEKKNRRQKH HRKSQRPSRT CQQFLKRCQL ASFTLPL
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	11.5 kDa
Expiration Date	6 months from date of receipt.