


Product Datasheet

Rat Eotaxin-2 protein (orb107708)

| | |
|----------------------------|---|
| Description | Recombinant Rat Eotaxin-2/CCL24. CCL24, also named MPIF-2, Eotaxin-2 and Ckbeta6, is a novel CC chem... |
| Conjugation | Unconjugated |
| Tested Applications | HPLC, SDS-PAGE |
| Target | Eotaxin-2 |
| Preservatives | Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4. |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only |
| Protein Sequence | MPTGSVTIPS SCCVTFISKK IPVNRVISYQ LANGSICPKA GVIFITKKGH KICTDPKLPWVQKHIKNLDA KRNQPSEGA ALGPKFVIQK LRGNSTKV |
| Purity | > 96% by SDS-PAGE and HPLC analyses |
| MW | 10.5 kDa |
| Expiration Date | 6 months from date of receipt. |



SDS-PAGE analysis of Rat Eotaxin-2 prote...