


## Product Datasheet

### FGF-4 protein (orb107712)

<b>Description</b>	FGF-4 is a heparin binding growth factor that is a member of the FGF family. Proteins of this family...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC, SDS-PAGE
<b>Target</b>	FGF-4
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of FGF-4 protein