


Product Datasheet

Murine FGF-9 protein (orb107714)

| | |
|----------------------------|---|
| Description | Fibroblast growth factor 9 (FGF-9) is a member of the fibroblast growth factor (FGF) family which is... |
| Conjugation | Unconjugated |
| Tested Applications | HPLC, SDS-PAGE |
| Target | FGF-9 |
| Preservatives | Lyophilized from a 0.2 µm filtered concentrated solution in 20mM Tris, 500mM NaCl, pH8.5. |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only |
| Protein Sequence | MPLGEVGSYF GVQDAVPFGN VPVLPVDSPV LLNDHLGQSE AGGLPRGPAVTDLDHLKGIL RRRQLYCRTG FHLEIFPNGT IQGTRKDHSR FGILEFISIAVGLVSIRGVD SGLYLG MNEK GELYGSEKLT QECVFREQFE ENWYNTYSSNLYKHVDTGRR YYVALNKDGT PREGTRTKRH QKFTHFLPRP VDPDKVPELYKDILSQS |
| Purity | > 95% by SDS-PAGE and HPLC analyses |
| MW | 23.3 kDa |
| Expiration Date | 6 months from date of receipt. |



SDS-PAGE analysis of Murine FGF-9 protei...