

Product Datasheet

Murine FGF-10 protein (orb107716)

Description	Recombinant Murine Keratinocyte Growth Factor-2/FGF-10. Fibroblast growth factor 10 belongs to the f...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-10
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, 600mM NaCl, pH7.4, 1mMmercaptoethanol
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QALGQDMVSQ EATNCSSSSS SFSSPSSAGR HVRSYNHLQG DVRWRRLLFSFTKYFLTIEKN GKVSQTKNED CPYSVLEITS VEIGVVAVKA INSNYLAMNKKGKLYGSKE FNNDCKLKER IEENGYNTYA SFNWQHNGRQ MYVALNGKGAPRRGQKTRRK NNTSAHFLPMT IQT
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	19.5 kDa
Expiration Date	6 months from date of receipt.