


Product Datasheet

Rat FGF-18 protein (orb107721)

Description	Recombinant Rat Fibroblast Growth Factor-18. Rat FGF-18 is encoded by the FGF18 gene in rats. By phy...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-18
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 500mM NaCl
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	EENVDFRIHV ENQTRARDDV SRKQLRLYQL YSRTSGKHIQ VLGRRISARG EDGDKYAQLLVETDTFGSQV RIKGKETEFY LCMNRKGKLV GKPDGTSKEC VFIEKVLNEN YTALMSAKYSGWYVGFTHKKG RPRKGPKTRE NQQDVHFMKR YPKGQTELQK PFKYTTVTKRSRRIRPHTPG
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	21.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat FGF-18 protein