


Product Datasheet

Rat FGF-21 protein (orb107723)

Description	Fibroblast growth factor 21 (FGF-21) encoded by the FGF-21 gene belongs to the large FGF family and ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	FGF-21
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, 500 mM NaCl, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	<p> AYPISDSSPL LQFGGQVRQR YLYTDDDQDT EAHLEIREDG TVVGTahrSPESLLELKALK PGVIQILGVK ASRFLCQQPD GTLYGSPHFD PEACSFRELLKDGYNVYQS EAHGLPLRLP QKDSQDPATR GPVRFLPMPG LPHEPQEQPGLPPEPPDVG SSDPLSMVEP LQGRSPSYAS </p>
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	19.7 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat FGF-21 protein